

Tissue Submission Memo Form

**FOR CLINIC TO FILL IN**

Requesting Date: Ordering Clinic Information/Stamp: Ordering Physician Name/Stamp:	Patient Information/Sticker Test Ordered: <input type="checkbox"/> APEX Tissue v.2 <input type="checkbox"/> COMPASS Tissue 1021 Genes <input type="checkbox"/> Idylla MSI Test
--	--

FOR PATHOLOGY LAB TO FILL IN

Dear Pathology Lab,

Please kindly prepare the FFPE slides according to the requirements below. Thank you.

- **1 H&E stained slide** with tumour region 5 - 25 mm² marked out, and tumour cellularity for **APEX & MSI ≥25%** and **COMPASS ≥20%**.
- **APEX/COMPASS: Matching 10 unstained sections** (Tumor area > 25 mm²) or **15 unstained sections** (Tumor area 5 - 25 mm²) of 5 µm thickness on uncoated/uncharged slides.
- **MSI:** Total tumor area 50 mm² [**Matching 2 unstained sections** of tumor area ≥ 25 mm² or max up to **Matching 5 unstained slides** of tumor area ≥ 10 mm²] of 5 µm thickness on uncoated/uncharged slides
- Histology report must accompany the slides.

To be completed by Pathology Lab

Specimen ID	FFPE Block ID										
Original specimen COLLECTION date	D	D	/	M	M	M	/	Y	Y	Y	Y
COLLECTION time	H	H	:	M	M	HR 24-hr time <input type="checkbox"/> No time info available					
Specimen source (e.g. Left lung)											
1 H&E stained slide	<input type="checkbox"/> Tumour region of interest (ROI) is 5-25 mm ² . <input type="checkbox"/> Tumour region of interest (ROI) is > 25 mm ² . Specify tumour cellularity in ROI: _____ %										
Unstained slides	No. of slides: ____ slides Thickness: ____ µm (≥5µm required)										
Slides SECTIONED date	D	D	/	M	M	M	/	Y	Y	Y	Y
Histology report	<input type="checkbox"/> Histology report attached										

Please send the slides [to the Ordering Clinic](#) once the slides are prepared.

Please call M Dx lab at +65-6950 3210 if you need any clarification.

M Diagnostics Lab. A subsidiary of MirXES Pte Ltd.
 30 Biopolis Street, #09-05/06 Matrix, Singapore 138671